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GB JC-362

JC-362 BODY MASSAGER

Thank you for purchasing this new Joycare body massager, ideal for slimming down, invigorating and massaging your body. This product is suitable for every part of your body, its oscillating and massaging heads penetrating deep down to help fight the effects of fluid retention.

For the product's correct usage, we recommend reading the following instructions carefully, and saving this operating manual for future use.

IMPORTANT SAFETY INSTRUCTIONS

Before connecting the appliance, check to make certain the voltage specified on the information nameplate corresponds to that of your power outlet. Before using the unit and its accessories, check to make certain no visible signs of damage are apparent. Always control the condition of the appliance and its accessories, and if any malfunctions or damage to the accessories are noted, immediately interrupt usage and contact an authorised servicing centre.

Always detach the power cord when the unit is not in use, before proceeding with its cleaning, when changing accessories, and whenever a malfunction or operating abnormality is noted. Always switch off the power switch before detaching the power cord. Always pull the plug and not the cord out its socket. Do not place the appliance where it can easily fall into the washbasin or bathtub. Never immerse the appliance in water or any other liquid.

Never use the appliance in the bath or shower, and if the unit should fall in water **DO NOT REACH FOR IT OR TRY TO RETRIEVE IT!**


Never use the appliance when sleeping or in case of drowsiness. Never use or place the unit underneath covers or pillows, as it may overheat and cause a fire and/or electrical shock. Keep the appliance out of the reach of children. Use the product only for the purposes described in this instruction manual, and do not use accessories other than those indicated in this manual. This product is designed for domestic use only, and is not intended for commercial use or for use in the medical sector. Always check the condition of the appliance, its accessories, power cord and plug, and interrupt usage immediately if any malfunctions or damage to the accessories are noted, contacting an authorised servicing centre. Always keep the appliance far from sources of heat, direct sunlight, humidity or sharp objects. Keep the power cord far from sources of heat and never coil it around the appliance. Do not insert any type of object into the slots on the unit. Do not block the air inlets. Never leave the unit unattended when in use. Do not operate where aerosol (spray) products are being used or in an environment in which oxygen is being administered by any other equipment. To prevent any possibility of accident or injury to children, never leave the unit unattended in their presence. Keep the appliance out of the reach of children.

SPECIAL INSTRUCTIONS FOR USING THE BODY MASSAGER

Treatments using this appliance must be pleasant and relaxing; should any pain or discomfort appear when using the unit, immediately interrupt the treatment and consult your physician. Do not use the unit on or around irritated or particularly sensitive skin. If you have any doubts regarding your health condition, always consult your physician before using the appliance. In the following cases, always consult your physician before using the appliance: pregnant women, postoperative convalescence, postoperative wounds not yet healed or inflamed; serious pathologies, cardiac diseases, bone fractures, sprains, for anyone affected by physical anomalies, and in the event of skin swelling up or blisters. It is advisable to always start the massage at a low setting and gradually increase the intensity. It is advisable to make use of the unit one hour after meals. Do not use the unit for more than 15 minutes on the same part of your body. Do not use the unit on your head. Never use the unit before going to sleep: the stimulating effect of the massage can result in a temporary state of insomnia. Do not use the appliance for more than 15 minutes consecutively; this is necessary to prevent any possible overheating which would reduce the unit's service life. In the event that treatments exceeding 15 minutes are necessary, it is advisable to switch off the unit, allow it to cool down and then resume its usage.

COMPONENTS AND FUNCTIONS

1. On/Off switch and intensity and massage speed control
2. Handle
3. Soft disc: for a more relaxing massage, protective head (mesh accessory) for massaging areas in which body hair can get caught up in the unit
4. Abrasive disc: provides exfoliating action for feet and heels
5. Disc with 8 rotating pearls: for the application of slimming creams or moisturising lotions
6. Flat disc: ideal for massaging the neck, shoulders and back
7. Wave disc: with 4 small waves, for deep massaging action

Note: The mesh accessory must be used with the flat and wave head attachments only, and never with the roller or abrasive attachments.

Recommendations for the use of the various accessories: arms and underarms (toning up); shoulders (alleviates tension and stress on shoulders and back); abdominals (slims down and carves abdominal muscles); hips (models and tones up hips and buttocks); legs (slims down legs and alleviates tension); calves (relaxes sore, tired muscles); feet (removes calluses and dry skin).

USING THE UNIT

This appliance is ideal for slimming down, invigorating, alleviating muscular tension and improving blood circulation, providing a sensation of relief and relaxation while massaging every part of your body. Also favours the application and absorption of slimming creams and moisturising lotions.

How to insert and remove the accessories: remove the mesh accessory by pulling, then place your fingers on the two hollow areas and pull out; to insert the accessory simply push.

After having inserted the desired accessory, based on the above descriptions, connect the unit to the electrical socket, making certain beforehand that the power switch is in the OFF position. Switch on the unit and set the desired massaging intensity by turning the roller. When the massage is completed, set the power switch in the OFF position and disconnect the power cord.

MAINTENANCE AND STORAGE

Always detach the power cord before proceeding with cleaning the unit. Do not use any type of solvent to clean the appliance. Clean with a soft, slightly dampened cloth, and never use detergents. Never immerse the appliance in water or any other liquid. Never use chemical detergents. Always store the unit in a dry, well ventilated location, away from heat sources, direct sunlight, or sharp objects. Never coil the power cord around the appliance, and do not use the power cord to hang up the unit.

TECHNICAL SPECIFICATIONS

Description: JC-362 body massager

Massage intensity setting: up to 2500 rotations per minute

Non-slip grip

5 interchangeable diskettes for different uses

Power specifications: 220-240W - 50/60Hz

Power: 25W max.

Weight: 1 Kg

Material: Abs

Class II appliance

Do not use this appliance in the vicinity of a bathtub, shower, washbasin or other recipients containing water.

These devices conform to all applicable European directives.

DISPOSAL

At the end of its life cycle, the device (including its removable parts and accessories) must not be disposed of together with household waste, but in compliance with European Directive 2002/96/EC. Since it must be handled separately from household waste, either bring it to a recycling and waste disposal centre for electrical and electronic appliances or give it back to the retailer on purchasing a new equivalent unit. Any infringement will be severely prosecuted.

Specifications and designs are based on the latest information available at the time of printing, and are subject to change without notice.

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IMPORTANT NOTICE

This appliance is not intended for use by people (including children) with reduced physical, sensory or mental capabilities, or lack of experience or knowledge, unless they are supervised of have been given instructions on the use of the product by a person responsible for their safety. Children should be supervised to ensure they do not play with the appliance. If the power cord is damaged, it must be replaced by the manufacturer or by its technical service or, in any case, by a person with a similar qualification in order to avoid risks. This appliance has heating elements, and thus, those who are not sensitive to heat should pay special attention when using it. Unplug the appliance after use if operating it in the bathroom as the vicinity of water is potentially dangerous even when the appliance is turned off.

The packaging (plastic bags, cardboard, polystyrene, etc.) are potentially dangerous and must be kept out of the reach of children. They must be disposed of according to current regulations.

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GB: WARRANTY CARD

Joycare S.p.A., a company with a sole shareholder, offers a period of 3 year warranty after the date of purchase. The warranty does not cover damage caused by accidental breakage, transportation, incorrect maintenance or cleaning, misuse, tampering or repairs carried out by non authorised third parties. The following is also excluded: damage caused by incorrect product installation, and damage caused by wear and tear, obviously including parts subject to wear and tear. In order to obtain the replacement/repair of products under warranty in cases other than those mentioned above, but which are instead caused by manufacturing defects, please contact your local importer/distributor. Goods that are not held to be in good condition because of the aforesaid causes must be accompanied by this certificate, duly filled-in, as well as a receipt for the goods (invoices and transport documents) and purchasing documentation.

Please keep the original packaging of the product.

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JOYCARE S.p.A. socio unico